


AIDIAN



Near-patient SARS-CoV-2 testing.

Qlife egoo.health system

Say hello to Egoo!

Near-patient SARS-CoV-2 testing.

The egoo.health system is a small and transportable, automated in vitro diagnostic device for near-patient testing of biological samples. The egoo.health system utilizes disposable SARS-CoV-2 test cassette with Aidian's proprietary Strand Invasion Based Amplification (SIBA®) technology.

Easy-access molecular testing of SARS-CoV-2

- Automated SARS-CoV-2 cassette test
- High quality results in 30 minutes
- Less than 2 minutes of hands-on-time
- Transportable device
- Random access testing of samples (with individually controlled units)
- Can be used both in laboratory or remote testing sites

Innovative technology and design

- Scandinavian design
- Uses Aidian's proprietary Strand Invasion Based Amplification (SIBA) technology
- Extremely small device, weight only 0.48 kg
- Tiny footprint which saves desk space, dimensions 55 x 100 x 110 mm
- Easy to use software: Egoo Clinical application
- No cable connection to PC or tablet computer (Utilizes Wi-Fi connection)



For more details, please consult instructions for use.

Contact info

Aidian Oy

aidian@aidian.eu

+358 10 309 3000

For orders: orders@aidian.eu

www.aidian.eu

Product	Prod. no	Order no
Egoo Q200 DK	E78852	154845
SARS-CoV-2 Capsule 10 Box DK, 10 pcs	E48364	154843
Egoo SARS-CoV-2 Lysis Buffer system, 50 pcs	E70464	154844

AIDIAN